

AVICHE

Automatic Digital Blood Pressure Monitor

JAPANESE CHIP TECHNOLOGY / PRECISE AND DURABLE

ES10

ES10




JAPANESE CHIP TECHNOLOGY

More than 850 hospitals and clinics are using in China market.




Model	ES10
Display	LCD display with backlight
Voice broadcast	English
Power source	4 "AA" batteries 1.5V/USB charging cable
Measurement range	Pressure: 0-280mmHg (0-37.3kPa) Pulse: 40 to 199 beats / min
Accuracy	Pressure: ± 3 mmHg (± 0.4 kPa) Pulse: $\pm 5\%$ of display reading
Pressure sensor	Semiconductor
Inflation	Automatic by electric pump
Air release	Automatic exhaust valve
Memory Functions	Up to 60 groups (two-person memories)
Dimensions (L x W x H)	135 x 110 x 80 mm
Net weight	279g (not including batteries)
Cuff circumference	22-36cm or 26-42cm(optional)
Automatic power off	1 minute
Work mode	Continuous


- ✓ **Large Backlight LCD Display:** Easy to read large display, the back-light makes it easy to read even in the dark.
- ✓ **Voice Broadcasting intelligent voice**
- ✓ **Core upgrading japanese technology:** Japanese chip technology precise and durable
- ✓ **One-key operation:** Press on the power button, this monitor automatically measures systolic, diastolic, pulse and display the readings on screen.
- ✓ **Double power source mode:** Double power source mode USB Cable or 4XAA(1.5V) Batteries. Automatic shutdown 60 seconds after measurement, saving 90 % of power.




Two User Mode




± 3 mmHg Accuracy




Professional Design



Large LCD Display



Heartbeat



Fully Automatic

• Packing List



Blood pressure monitor




BP cuff



Usb Power Cord



Storage bag



User manual