


Automatic Digital ES12 Blood Pressure Monitor



Accuracy



Two User Mode



Backlit LCD Display



Irregular Heartbeat Indicator



Voice Broadcast



Fully Automatic





LCD Backlit Display: Easy to read display, the back-light makes it easy to read even in the dark.

Heart rate test + Core upgrading aviation technology: In the measurement of blood pressure at the same time on the heart rate is also tested, a machine more use.

One-key operation: Press on the power button, this monitor automatically measures systolic, diastolic, pulse and display the readings on screen.

Record Two People data: More Saving Money
2 x 100 Large-capacity memory values for easy comparison

Double power source mode: Double power source mode USB Cable or 4XAA(1.5V) Batteries. Automatic shut-down 60 seconds after measurement, saving 90 % of power.

Model	ES12
Display	English
Measurement Method	Oscillometric method
Power source	DC6V (4 “AAA” batteries 1.5V)
Pressure	0 to 280 mmHg (0kPa–37.5kPa)
Pulse	40 to 199 beats/min
Accuracy	Pressure: ± 3 mmHg (± 0.4 kPa) Pulse: $\pm 5\%$ of display reading
Inflation	Fuzzy-logic controlled by electric pump
Dimensions (L x W x H)	133 x 116 x 81 mm
Deflation	Automatic pressure release valve
Memory	Up to 100 per user
Time	Year, month, hour, minute
Operating Temperature/Humidity	5° C to 40° C/15% to 80%RH
Storage Temperature/Humidity	-20° C to 55° C/15% to 90%RH

PACKING LIST

