

AVICHE

Medical **5L** Oxygen Concentrator

It is a new generation oxygen concentrator for homecare and medical usage, which use our patented technology and have the ability to adjust the flow rate at 0.1L/min and moreover, there is a large display which can real-timely display oxygen concentration and oxygen flow. All the designs are customized to meet precise operation requirement and enhance the customer's experience.



Countinous
Work Mode



93% \pm 3% oxygen
concentration



Oxygen &
Nebulization



Precise
adjustment



Medical Oxygen Concentrator

Take care of family health with oxygen



1 Precise adjustment

0.1L
Mechanical adjustment
Oxygen flow
real-time display


Uses patented technology to adjust the oxygen flow at the rate of 0.1l/min, users can rotate the control button to adjust the oxygen flow.



2 High-performance compressor

Continuous oxygen
5 years Service life

Strong power, stable operation without pressure in high altitude areas, and the concentration is still maintained above medical standards



3 Stable Oxygen Supply

5L/min
Oxygen Output
93% ± 3%
Oxygen Concentration

High-density molecular sieve, Improved service life and more stable oxygen concentration



4 Nebulization

The oxygen and nebulization can be switched at any time to meet the choices of more families.

* Nebulization can help treat respiratory diseases such as cough, pneumonia, bronchitis, etc., non-toxic and side effects are safer



Oxygen flow	1.0-7.0L/min
Oxygen concentration(5L)	93% ± 3%
Nebulization	✓
Input power	220V 50Hz / 110V 60Hz

Weight	20.5kg
Dimension	W30cm x L37.5cm x H59.8cm
Wattage	<400VA
Noise	<46dB