

Medical **A15+A25** Oxygen Concentrator

A15 and A25 is a new generation oxygen concentrator for homecare and medical usage, which use our patented technology and have the ability to adjust the flow rate at 0.1L/min and moreover, there is a large display which can real-timely display oxygen concentration and oxygen flow. All the designs are customized to meet precise operation requirement and enhance the customer's experience.



Medical Oxygen Concentrator

OEM/ODM

Large LCD Display

Large touch screen make it easy to use.



Oxygen concentration: >90%

Adopted advanced aviation technology and use double molecular sieve to make oxygen, high oxygen concentration, steady running and durable.



Precise Oxygen Flow Adjustment

Patented technology to adjust the oxygen flow at the rate of 0.1 Liter/min. Users can just press the control button to adjust the oxygen flow.



● Application



Can be used in pet oxygen cabins

10 million pcs/cm³ negative ions

Equipped negative IONS generator, release 10 million pcs/cm³ negative ions, purify inlet air and keep oxygen more fresh.



6M Distance Remote Controller

Users can remotely control the oxygen concentrator even lying on the bed or sitting on a chair. It is more easy and conformable to use.



● Product Parameters



A15

Oxygen flow : 1.0-3.0L/min
 Oxygen concentration(1L): >92%
 Input power:220V/110V
 50Hz/60Hz
 Wattage: 130W
 Weight.: 9.1kgs
 Maximum noise level:<50dB
 Size:302 x 278x 375 mm



A25

Oxygen flow : 1.0-5.0L/min
 Oxygen concentration(2L): >92%
 Input power:220V/110V
 50Hz/60Hz
 Wattage: 160W
 Weight: 9.1 kgs
 Maximum noise level:<50dB
 Size:302 x 278 x 375mm